

Harnessing data to transform the power of molecular diagnostics and epidemiological modelling

Jesus Rodriguez Manzano, PhD


Lecturer in Antimicrobial Resistance and Infectious Disease at Imperial College London
Deputy Director of the Centre for Antimicrobial Optimisation (CAMO)

Analysis and prediction of carbapenem resistance

Outbreak Surveillance



Spatial-temporal inferences



Ashleigh Myall

Computational
Researcher

Predicting Outbreaks

$$f(\text{hospital}) = \text{virus particles}$$

Acquisition/Res Prediction

$$f(\text{person}) = \text{person with virus}$$

Research topics:


Graph-based machine learning, Spatial-temporal data, Decision support systems & antimicrobial-resistance.

Funders:




Data-driven multiplexing: concept

Maintain traditional qPCR workflow and instrumentation while streamlining data analysis and management




Data-driven multiplexing: how does it work?



Accuracy: **99.1%**
(20% increase)



True label	bla_{NDM}	4558	5	26
	bla_{VIM}	2	5623	57
	bla_{KPC}	9	49	5859
		bla_{NDM}	bla_{VIM}	bla_{KPC}
		Predicted label		



Amplification and melting curve analysis



AMPLIFICATION CURVE ANALYSIS (ACA)

Kinetics information



MELTING CURVE ANALYSIS (MCA)

Thermodynamics information



Accuracy

ACA

MCA


82.31%

89.34%

Amplification and melting curve analysis

Classification accuracy:

99.28%




Data-driven multiplexing: CPE as case study

Target	N	TP	TN ^a	FP	FN	SEN	SPE	Accuracy (CI)
<i>bla</i> _{IMP}	45	45	32	0	0	100.0%	100.0%	100.0% (95.32 to 100.00%)
<i>bla</i> _{KPC}	9	9	32	0	0	100.0%	100.0%	100.0% (91.40 to 100.00%)
<i>bla</i> _{NDM}	74	73	32	1 ^b	0	100.0%	96.97%	99.06% (94.86% to 99.98%)
<i>bla</i> _{OXA-48}	84	84	32	0	0	100.0%	100.0%	100.0% (96.87 to 100.00%)
<i>bla</i> _{VIM}	8	8	32	0	0	100.0%	100.0%	100.0% (91.19 to 100.00%)
<i>bla</i> _{OXA-48} & <i>bla</i> _{NDM}	1	1	32	0	0	100.0%	100.0%	100.0% (97.24 to 100.00)
Total	221	220	32	1	0	100.0%	96.97%	99.60% (97.82 to 99.99%)

Abbreviations - N: number of samples; TP: True Positive TN: True Negative; FP: False Positive, FN: False Negative, SEN: Sensitivity, SPE: Specificity, CI: Confidence Interval.

^aA total 32 negatives samples are considered across all the groups for sensitivity, specificity and accuracy calculation

- 7-plex for respiratory pathogens using TaqMan probes & LAMP (*Sensors & Diagnostics*. 2022)
- 5-plex for CPE using Sybr green (*Anal Chem*. 2019 & *Front Mol Biosci*. 2021)
- 9-plex for mobilized colistin resistance genes using Sybr green (*Anal Chem*. 2020)



camo
centre for
antimicrobial
optimisation



NIHR | National Institute
for Health Research



UK
RI UK Research
and Innovation

Rosetrees Trust
Supporting the best in medical research

Community
Jameel

